

# ELECTRONIC PNEUMATIC TOURNIQUETS

EASY PUMP - EASY PUMP DUAL

Marking	Class
CE	Ila

## EASY PUMP DUAL

Ref.G10903



## EASY PUMP

Ref.G10803



- Touch screen interface : 120mm x 90mm
- 1 or 2 regulated pressure circuits
- Setting range : 0 - 600 mmHg
- Controlled by microprocessor
- Battery
- Compact


**TRACEABILITY**



Alphanumeric keyboard

**SECURITY**

Alarms	Warning
Pressure	Audible and visual
Time	Audible and visual
Battery	Audible and visual



**EASE OF MAINTENANCE :**

- All fuses are externally accessible
- Removable battery via an access hatch



Dimensions (in mm)	Height	250
	Width	300
	Depth	190
Weight (in Kg)	5.8	

Mobile stand on casters  
**Ref.A10701**



Printer  
**Ref.G13200**



# ELECTROMECHANICAL PNEUMATIC TOURNIQUETS

## LITTLE PUMP - LITTLE PUMP DUAL

Marking	Class
CE	Ila



### LITTLE PUMP

Ref.G10705

#### Easy to use

- 1 regulated pressure circuit
- Controlled by microprocessor, a guarantee of security and accuracy
- Setting pressure and timer
- Setting range : 0 - 600 mmHg, with an accuracy of  $\pm 5$ mmHg
- Timing and pressure alarms
- **Optional** : Mobile stand on casters **Ref.A10701**

Dimensions (in mm)	Height	238
	Width	167
	Depth	160
Weight (in Kg)	2.2	

Clip for hanging on infusion stand **Réf.A10760**



**LITTLE PUMP DUAL**  
**Ref.G10706**

- 2 regulated pressure circuits
- Controlled by microprocessor, a guarantee of security and accuracy
- Setting pressure and timer
- Setting range : 0 - 600 mmHg, with an accuracy of +- 5mmHg
- Timing and pressure alarms

Mobile stand on casters  
**Ref.A10701**



Dimensions (in mm)	Height	250
	Width	300
	Depth	190
Weight (in Kg)	4.2	



# MANUAL PNEUMATIC TOURNIQUET

MONOBLOCK - Decontaminable

Marking	Class
CE	I



Reference	Designation	Dimensions incm		
		Cuff width	Cuff length	Overall length
GMBM01	New born	3.5	29	38
GMBM02	Child	4.5	35	48
GMBM03	Small adult	6	46	61
GMBM04	Adult	8	46	67
GMBM05	Conical arm cuff Adult L	10	47	70
GMCM04	Conical lower limb cuff Adult	10	62	76
GMCM05	Conical lower limb cuff Adult L	10	76	90
GMCM06	Conical lower limb cuff Adult XL	10	82	100
GMCM07	Conical lower limb cuff Adult XXL	10	107	122
KGMBM1	1 arm cuff Child, adult and 1 lower limb cuff adult L	X	X	X
KGMBM2	1 arm cuff Adult and 1 lower limb cuff adult L	X	X	X

## STERILIZABLE

Silicone bladder

Ref.	Designation
G20154	Manual tourniquet arm cuff/lower limb cuff New born
G20153	Manual tourniquet arm cuff Child
G20155	Manual tourniquet arm cuff Small Adult
G20152	Manual tourniquet arm cuff Adult
G20150	Manual tourniquet lower limb cuff Adult
G20151	Manual tourniquet lower limb cuff Adult L
G20156	Manual tourniquet lower limb cuff Adult XL



## EMERGENCY PNEUMATIC TOURNIQUET

- Single patient use
- Safety with pressure indicator
- Maximum pressure 300 mmHg (safety valve)
- With protection for the skin

Reference	Designation	Dimensions in cm		
		Cuff width	Cuff length	Overall length
GMBU01	New born	3.5	29	48
GMBU02	Child	4.5	35	48
GMBU04	Adult	8	46	67
GMCU05	Conical lower limb cuff Adult L	10	47	70
KGMBU1	Kit with : 1 x Child, 1 x Adult, 1 x Conical lower limb cuff Adult L	X	X	X



# EASY TEST

Ref.ET903

Simply perform the maintenance of all our tourniquet :

- Touch screen interface : 120mm x 90mm
- Controlled by microprocessor
- Battery
- Compact

Marking	Class
CE	Ila

**New product  
2019**



### Compatible devices :

- G10903 - EASY PUMP DUAL
- G10803 - EASY PUMP
- G10705 - LITTLE PUMP
- G10706 - LITTLE PUMP DUAL
- MODULO FUSE

Control and verification	YES	NO
Sensor calibration	X	
Pressure in the circuits	X	
Sealing of pressure circuits	X	
Managing device pressure	X	
Human - machine interface	X	
Operating file	X	
Supply voltage	X	
Date and hour	X	